



NHLS Stool Parasites PT Scheme _____

PTS LAB NO: _____

LABORATORY NAME: _____

CHALLENGE:	ANSWER CODE/S:	CLINICAL HISTORY:	INSTRUCTIONS:
PS0 /21 Stool/urine concentrate			Vortex/mix the specimen well. Make a few wet preparations and examine the slide for parasites using 10x and 40x objectives.
PS0 /21 Stool/urine concentrate			Vortex/mix the specimen well. Make a few wet preparations and examine the slide for parasites using 10x and 40x objectives.
PS0 /21 Stool/urine concentrate			Vortex/mix the specimen well. Make a few wet preparations and examine the slide for parasites using 10x and 40x objectives.
PS0 /21- Stool/urine concentrate			Vortex/mix the specimen well. Make a few wet preparations and examine the slide for parasites using 10x and 40x objectives.
PS /21 Stool smear			You are provided with a fixed stool smear; stain with an appropriate stain and examine the slide for parasites

	Stain used:		using the appropriate objective. Please do not return your slide.
--	-------------	--	--

Email: parapts@nhls.ac.za or Fax: 086 225 2460

See closing date on the schedule on **www.nhls.ac.za**