

NMAL VIROLOGY LABORATORY –
SPECIFICATIONS – PIPETTES STAND WITH UV LIGHT

1. Pipette stand to fit all pipette brands from 3 - 6 single-channel pipettes.
2. Rotate freely and easy access to all pipette positions.
3. To have germicidal UV Lamp that disinfect and decontaminate pipette for about 15 – 30 minutes to reduce contamination.
4. UV: 220 - 300nm UV radiation to eliminate surface contamination with 450 micro watts.
5. Button to activate the lamp.
6. UV resistant polycarbonate.
7. Aluminium coated, polyester film.
8. Dimensions : 15-18cm diameter, 30 - 40cm tall & rotate freely
9. Electrical : 100V to 240V, 50/60 Hz, 1A
10. Easy to dissemble.
11. Warranty period – minimum 1 year.
12. Operators manual.

QA OFFICER: X. MADIKANE



24/06/2022

ENDUSER/LAB MANAGER: H. VILAKAZI



24/06/2022