

EQUIPMENT**Cytospin Centrifuge**

Business Unit	
Laboratory	
RFQ0008315	
Specification	Requirement
Foot print	Benchtop
Sealed head	Removable and autoclavable
Design	Window to enable viewing of sealed head
	Mechanical & electrical components protected from accidental spills
	Easy cleaning and decontamination
Control Panel	LCD display touch pad
Speed	Adjustable 200-2000 rpm
RPM/RCF mode	RPM
Acceleration	Variable acceleration
Deceleration	Variable deceleration
Programmability	Pre selection of speed and time
Specimen chambers	Minimum 12 samples at one time
	Compatible with filter cards with 2 circles
Cytoclips (metal)	Double
Cytofunnels (plastic)	Double
Run Time	Variable
Decibels	≤70 dBA
Ambient temperature	10 - 35°C
Internal temperature	Not applicable
Electrical	Compliant to South Africa
Delivery time within 6 weeks	Yes
In use in NHLS	Yes
Additional Requirements	
Manual opening when power failure	Yes
Electronic key pad	Yes
Door safety lock	Yes
Fast stop option	Yes
Ergonomics	Allows for one-handed opening & closing
Operators' manual	English
Installation and commissioning	By Supplier
Laboratory Manager	
Signature	
Date	
Quality Assurance	
Signature	
Date	