

Cecilia Makiwane NHLS Lab

Specifications for phlebotomy training arm

Purpose for use: Phlebotomy Practice Kit, Venipuncture Practice Arm, IV Practice Arm Kit for Nurses Practice & Training, IV Starter Kit, for Educational Training Teaching.

Packaging dimension: 23.6 x 6.9 x 5.7 inches (59.9 x 17.5 x 14.6 cm); 4.15 Pounds (1.87 kg)

Features and benefits:

Adult size practice arm

Realistic practice for IV injections

Withstand hundreds of repeated punctures without leakage

Used in social training institutions, clinics, hospitals.

Made of premium quality

Feels real to the touch and has anatomically correct cephalic or basilic veins with multiple puncture sites on the dorsal hand, lateral pulse, and the cubital fossa area on the ventral fore arm.

Automatically reseal after every stick

Adjusted heart rate and pulse strength simulating a heart rate from 10 - 150BPM

Ease to assemble

Cover multipurpose venous training arm, advanced venous training arm and arterial and venous IV arm with pulse and pump

Additional accessories

Phlebotomy arm user guide for training

2 x Adjustable phlebotomy arm rest: moveable

Laboratory Manager *BABONE* date *21/06/23* signature..... 